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MED<sup>9</sup>EL

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hearLIFE

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- Anderson I, Baumgartner WD, Böheim K, Nahler A, Arnolder C, D'Haese P (2006) Telephone use: What benefit do people with a cochlear implant receive? *International Journal of Audiology*, 45: 446-53.
- Carmel E, Kronenberg J, Wolf M, Migirov L (2011): Telephone use among cochlear implanted children. *Acta Otolaryngol*; 131(2):156-9.
- Castro A, Lassaletta L, Bastarrica M, Prim MP, De Sarriá, Gavilán J (2006) Comparison of different mobile telephones in COMBI 40+ users, *Acta Otolaryngol*, 126:714-7
- Dillon H, James A, Ginis J (1997) Client oriented scale of improvement (COSI) and it's relationship to several other measures of benefit and satisfaction provided by hearing aids, *Journal of American Academy of Audiology*, 8(1): 27-43
- Erber NP (1985) Telephone communication and hearing impairment, Taylor and Francis: London
- Giles E (2005) Telephone training with a cochlear implant in Pedley K, Giles E, Hogan A (Eds) *Adult cochlear implant rehabilitation*, Whurr: London
- Ricability, Stay in touch: A guide to telephones and services for older and disabled people, [www.ricability.org.uk](http://www.ricability.org.uk)
- Spitzer JB, Leder SB, Giolas TG (1993) Rehabilitation of late-deafened adults: Modular program manual, Mosby: St Louis
- Tait M, Nikolopolous TP, Archbold S, O'Donoghue GM (2001) Use of the telephone in prelingually deaf children with a multichannel cochlear implant, *Otology and Neurotology*, 22: 47-52
- Wayner DS (1998) *Hear what you've been missing: How to cope with hearing loss*, John Wiley Publishing: New York
- Wayner DS, Abrahamson JE (1998) *Learning to hear again with a cochlear implant; Clinician Manual; A Personal Journal; User's Guide*, Hear Again Inc: New York

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MED-EL UK Head Office  
office@medel.co.uk

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MED-EL Colombia S.A.S.  
Office-Colombia@medel.com

MED-EL Indonesia  
office@id.medel.com

VIBRANT MED-EL  
Hearing Technology France  
office@fr.medel.com

MED-EL Mexico  
Office-Mexico@medel.com

MED-EL Korea  
office@kr.medel.com

MED-EL BE  
office@be.medel.com

MED-EL Middle East FZE  
office@ae.medel.com

MED-EL Vietnam  
office@vn.medel.com

MED-EL GmbH Sucursal España  
office@es.medel.com

MED-EL India Private Ltd  
implants@medel.in

MED-EL Japan Co., Ltd  
info@medel.co.jp

MED-EL GmbH Sucursal em Portugal  
office@pt.medel.com  
MED-EL UK Ltd

MED-EL Hong Kong  
Asia Pacific Headquarters  
office@hk.medel.com

MED-EL Liaison Office Australasia  
office@au.medel.com

MED-EL Philippines HQ  
office@ph.medel.com

MED-EL Medical Electronics  
Fürstenweg 77a | 6020 Innsbruck, Austria | office@medel.com

medel.com